

Cytology



Cell sampling solutions...

US Endoscopy's line of cytology devices addresses a variety of needs when sampling from the upper and lower GI tract.

Infinity[®] cytology device

The premium choice, the Infinity[®] cytology device is designed for collecting substantial and quality samples.

- features a unique bristle configuration
 - stiffer bristles to help create a defect in the tissue
 - softer bristles to capture abraded material
- designed with spaces between the bristles to help pack cells into the catheter

Cytology brushes

These brushes are designed to collect cell samples in the bronchus as well as the GI tract.

- features a rotatable handle to facilitate an abundant sample

For more information on cytology devices, contact US Endoscopy today.

| product number | description | sheath diameter (mm) | length (cm) | units/box |
|----------------|---------------------------------------|----------------------|-------------|-----------|
| 00711401 | bronchoscope | 1.8 | 120 | 20 |
| 00711402 | gastroscope | 1.8 | 160 | 20 |
| 00711403 | colonoscope | 2.4 | 230 | 20 |
| 00711405 | enteroscope | 2.4 | 350 | 10 |
| 00711499 | Infinity [®] cytology device | 2.4 | 230 | 5 |

Distributed by:

US endoscopy

listening...and delivering solutions[®]

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com