

Oracle™
EUS balloon



listening... and delivering solutions®

A durable EUS balloon offering superior tear resistance...

The Oracle™ EUS balloon is manufactured not to contain latex and allows for easy application.

The Oracle™ EUS balloon...

- minimizes balloon tearing during application to the endoscope
- eliminates the risk of exposing patients and staff to latex
- an elastic band design allows for easy application
- elastic material offers a secure grip to ensure the balloon stays on the endoscope throughout the procedure



radial array

One balloon for both Radial and Curvilinear endoscopes.



curvilinear array

For more information on the Oracle™ EUS balloon, contact US Endoscopy today.

product number	description	sterile	units/box
00711891	Olympus curvilinear/ radial array	yes	20

Distributed by:

US endoscopy

*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

760868 Rev. A

Olympus® is a registered trademark of Olympus Optical, Inc.
US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy