

**Cheirón**<sup>®</sup>   
...we breathe for you.

EU EACH PUMP ✓  
**TESTED**  
SER. No 8888 SINCE 1994



# VICTORIA

ICONIC SUCTION UNIT

- **EXTREMELY SILENT RUN**
- MAINTENANCE-FREE GERMAN AGGREGATE
- AGGREGATE BODY IS MADE OF ALUMINIUM ALLOY
- RELIABLE AND VIBRATION-FREE RUN
- AGGREGATE IS COOLED BY EXTRA FAN
- DESIGNED FOR CONTINUOUS RUN
- LOW MAINTENANCE REQUIREMENTS
- **6 YEAR WARRANTY ON THE AGGREGATE**
- INNOVATIVE EUROPEAN DESIGN

[www.cheiron.eu](http://www.cheiron.eu)

- **EXTREMELY SILENT RUN**
- MAINTENANCE-FREE GERMAN AGGREGATE
- AGGREGATE BODY IS MADE OF ALUMINIUM ALLOY
- RELIABLE AND VIBRATION-FREE RUN
- AGGREGATE IS COOLED BY EXTRA FAN
- DESIGNED FOR CONTINUOUS RUN
- LOW MAINTENANCE REQUIREMENTS
- **6 YEAR WARRANTY ON THE AGGREGATE**
- INNOVATIVE EUROPEAN DESIGN



Overflow protection

Vacuum gauge in kPa, mmHg and Bar, colour coded

Illuminated on/off switch meets EU standard EN ISO 10079-1

Ergonomic vacuum regulator

Casing is made of very sturdy, unbreakable, fire retardant material with UV protection

Change-over valve shift the vacuum from one jar to the other with a single move – optional

Jar holder is equipped with suction circuit holder

The stainless eurorail ensures maximum flexibility. You can attach up to four jars of any size and/or any other accessories needed

## MAKES LESS NOISE THAN YOUR PC!

A SPECIAL SYSTEM OF NOISE AND VIBRATION DAMPENING FOR AN EXTREMELY QUIET AND SMOOTH OPERATION WITH AN INCREDIBLY SILENT RUN

Non-marking casters with safety brakes.  
Antistatic casters - optional

## VICTORIA VERSA

The most popular high performance model of Victoria especially suitable for wards, operating rooms, recovery rooms, ICUs, obstetrics and endoscopic applications, etc.

## VICTORIA LIPOS

High performance solution designed particularly for liposuction and operating rooms.

## VICTORIA PORTABLE

Quiet and effective suction unit for medical practices, wards and endoscopic applications, as well suitable for fast suction in any situation.

## VICTORIA THORAX

A suction unit specially designed for common use as well as for thoracic and abdominal surgery drainage and intensive care. The system can be used for long-term up to 30 days continuous drainage.

Thanks to a unique patented solution Victoria Thorax enables full suction output even during drainage.



## VICTORIA PORTABLE WITH EURORAIL

Advantage of the unit is that instead of regular holder, unit is equipped with eurorail 25 x 10 mm this allows you to use your holders and jars.



## TECHNICAL SPECIFICATIONS

TYPE	PORTABLE	VERSA	THORAX	LIPOS
Suction jar 2 l with press on cover	•	•		•
Suction jar 2 l with screw cover			• 2 pcs	
Suction jar holder for 2 l jar for the side of the unit	•			
Suction jar holder for 2 l jar for eurorail		•	• 2 pcs	•
PVC suction circuit with stop valve	•	•	•	•
Mobile stand		•	•	•
Limiting water valve			•	
Suction capacity nominal +/- 5 %	43 l/min	43 l/min	43 l/min	50 l/min
Suction capacity effective +/- 5 %	38 l/min	38 l/min	38 l/min	45 l/min
Vacuum	0-90 kPa (0-680 mmHg)	0-90 kPa (0-680 mmHg)	0-90 kPa (0-680 mmHg)	0-93 kPa (0-700 mmHg)
Drainage			0 - 40 cmH <sub>2</sub> O	
Input Power	110 VA	110 VA	110 VA	110 VA
Input Voltage	230V /50-60 Hz and models for 120V ± 10% / 60Hz are also available			
Dimensions (w x h x l) in millimeters	440 x 300 x 310	480 x 900 x 460	480 x 900 x 460	480 x 900 x 460
Dimensions (w x h x l) in inches	20,1 x 11,8 x 13,4	18,1 x 18,9 x 39	18,1 x 18,9 x 39	18,1 x 18,9 x 39
Weight (without accessories)	12 kg (26,5 lb)	21 kg (46,3 lb)	24 kg (52,9 lb)	21 kg (46,3 lb)

# OPTIONAL ACCESSORIES

Vacuum gauge in kPa, mmHg and Bar, colour coded



Bottle against overflow Victoria

Order No: 0-52-0002



Support ring for 2 l jar

Order No: 11-5141



Clamp 25x5

Order No: 11-5122



New eurorail for placing a large number of optional accessories

Order No: 11-5104



Disposable secreta bag MONOKIT

Order No: 114-087-001



0,5 – 4 l polycarbonate and polysulfone suction jars, autoclavable

Order No: Price list



MSF microbiological suction filter (20pcs/box)

Order No: 271-022-001



Easy and safe handling



Change-over valve

Order No: 11-5107



IEC power outlet. It is possible to order unit with IEC power outlet



Air-operated on/off footswitch

Order No: 114-047-000



Shelf

Order No: 11-5116



Storage basket

Order No: 11-5105



Holders allow you to place extra jars



Suction circuit holder

Order No: 11-5110

