

134 °C
273 °F

Biliary Dilation Catheter

endoscopic usage


- radiopaque teflon tube
- accepts guide wire .035 inch



| length | catheter Ø | catheter tip | colour mark | Order No. |
|--------|------------|--------------|-------------|-----------|
| 180 cm | Fr. 6.0 | Fr. 4.0 | violet | 200 00064 |
| 180 cm | Fr. 7.0 | Fr. 5.0 | yellow | 200 00075 |
| 180 cm | Fr. 8.5 | Fr. 5.0 | orange | 200 00085 |
| 180 cm | Fr. 9.0 | Fr. 5.0 | white | 200 00095 |
| 180 cm | Fr. 10.0 | Fr. 5.0 | red | 200 00105 |
| 180 cm | Fr. 10.0 | Fr. 7.0 | red | 200 00107 |
| 180 cm | Fr. 12.0 | Fr. 7.0 | blue | 200 00127 |

recommended equipment

Appropriate guide wires can be found in chapter 06.



reusable

Biliary Dilation Catheters