

efficiency

BioShield® irrigator...

Flush blood, exudate and other debris out of your way at anytime throughout a procedure*- even with a device in place

How important is a clear view?

The BioShield® irrigator affords immediate, direct, intraprocedural irrigation. Simply attach a 60cc syringe to the irrigation line, or use extension tubing to connect with your water pump for hands-free, foot pedal control.



1133

0071**1133**
Olympus/Fujinon

0071**1137**
Pentax (not pictured)

Add a whole new level of efficiency to your endoscopy procedures with the BioShield® irrigator

*immediate, direct irrigation
with device in channel*

*In order to achieve moderate irrigation pressure and minimize leakage, there must be sufficient room within the biopsy channel for fluid to pass when there is a device in place

convenience

product information

	part number	description	length (cm)	qty/box
Olympus & Fujinon (G5 series and newer) compatible				
	00711124*	biopsy valve (individually packaged)	—	200
	00711125	biopsy valve (20 per pack/10 packages)	—	200
	00711128	biopsy valve sterile (individually packaged)	—	100
	00711135	biopsy valve (individually packaged)	—	—
	00711126	biopsy valve (individually packaged)	—	100
	00711129	biopsy valve (20 per pack/10 packages)	—	200
	00711133	irrigator	30.5	50
Pentax compatible				
	00711127	biopsy valve (individually packaged)	—	100
	00711136	biopsy valve (individually packaged)	—	100
	00711137	irrigator	30.5	50
universal irrigation accessories				
	00711131	universal irrigating adaptor	2.54 (adaptor)	50
	00711134	universal water pump extension tubing	180	50

*Not currently for sale in the United States