

bite blocks



A broad offering to fit all patient needs...

Our extensive line of bite blocks offers a variety of options in sizes, designs, and models to

Standard

The sturdy, yet forgiving material of the standard bite blocks protects delicate instruments and provides comfort for patients.



maxi bite block / 60FR

- available with a velcro, latex-free or woven strap, or without a strap
- with dental retention rim



standard bite block / 48FR

- available with velcro and latex-free straps or without a strap
- with dental retention rim



intermediate bite block / 48FR

- latex-free strap
- without dental retention rim

product number	description	opening size	dental retention rim	units/box
00712801	standard - latex-free strap	16mm / 48 FR	yes	100
00712802	standard - no strap	16mm / 48 FR	yes	100
00712803	standard - velcro strap	16mm / 48 FR	yes	100
00712804	maxi - latex-free strap	20mm / 60 FR	yes	100
00712805	maxi - no strap	20mm / 60 FR	yes	100
00712806	maxi - velcro strap	20mm / 60 FR	yes	100
00712809	intermediate - latex-free strap	16mm / 48 FR	no	100
00712812	maxi - woven strap	20mm / 60 FR	yes	100

New

ensure the best for the patient and type of procedure.

Specialty

From a pediatric bite block to a purpose-built O₂ delivery bite block, our specialty bite blocks are designed to fit the unique needs of patients.



pediatric bite block / 38FR

- latex-free strap
- without dental retention rim



Vasto™ bite block / 72FR

- Featuring larger side ports for scope repositioning and suction access
- Latex-free strap
- With dental retention rim



Respa™ O₂ delivery bite block / 60 FR **New**

- Can be used to house the nasal cannula during procedures to easily deliver O₂ to the patient orally
- Latex-free strap
- With dental retention rim

For more information, contact US Endoscopy today.

product number	description	opening size	dental retention rim	units/box
00712808	pediatric - latex-free strap	12.6mm / 38 FR	no	25
00712810	Vasto™ bite block – latex-free strap	24mm / 72 FR	yes	100
00712811	Respa™ O ₂ delivery bite block	20mm / 60 FR	yes	100

New

Velcro® is a registered trademark of Velcro Industries B.V.

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy