

Ø 1.8  
mm

134 °C  
273 °F

single / double sided

- for cleaning of endoscope working channels



brush Ø 3.0 mm

length	single sided, incl. handle	double sided
120 cm	41 181-A	BS 41 181-A
160 cm	41 181-B	BS 41 181-B
180 cm	41 181-C	BS 41 181-C
280 cm	41 181-M	

Ø 2.2  
mm

134 °C  
273 °F

single / double sided

- for cleaning of endoscope working channels



brush Ø 5.0 mm

length	single sided, incl. handle	double sided
120 cm	41 221-A	BS 41 221-A
160 cm	41 221-B	BS 41 221-B
180 cm	41 221-C	BS 41 221-C
230 cm	41 221-G	BS 41 221-G
280 cm	41 221-M	

single / double sided

- for cleaning of endoscope working channels



brush Ø 7.0 mm

length	single sided, incl. handle	double sided
120 cm	41 401-A	BS 41 401-A
160 cm	41 401-B	BS 41 401-B
180 cm	41 401-C	BS 41 401-C
230 cm	41 401-G	BS 41 401-G