

ENDOTRANS **Motorised**

Endoscopy Procedure Cart



MEDICAL
INNOVATIONS
GROUP®

EndoTrans Features



Keyboard drawer
W: 410mm x D: 245mm



2 rear doors with 'push to open' handles



Charging point (included with motorised option)



Easy to use forward and reverse thumb control



Lockable drawer with secure combination lock and key



4 x twin castor wheels (2 braked)

EndoTrans is a highly versatile procedure cart and transportation solution supplied as standard with a large and easy to clean work surface, one lockable drawer, 2 equipment shelves and large 6" castors (2 lockable). Available in **motorised** and **non-motorised** versions, EndoTrans can be customised with optional accessories listed on page 5 & 6.

EndoTrans is compatible with most Olympus, Pentax and Fuji equipment. (Complete list of compatible equipment available on request).

Medical Innovations Group supplies custom procedure carts for many endoscopy companies around the world.

Product Ref	Product Description	Weight
104468	EndoTrans Procedure Cart - non-motorised	95kg
104469	EndoTrans Procedure Cart - motorised	102kg

Customised Options

EndoTrans is fully compatible with our MIG 30 Monitor cart which attaches securely using the pre-fitted 'quick connect and release' docking station as shown below. The MIG 30 cart includes a swivel & tilt LCD monitor mount and 2 large, easy grip pushing handles. Procedure cart and monitor cart move together when connected.



The MIG 30 Monitor Cart connects to EndoTrans via the easy-to-use docking station at the back of the unit



MIG 30 Monitor Cart

Product Ref	Product Description	Weight
104470	MIG 30 Monitor Cart includes LCD swivel & tilt (max weight 20kg), two pushing handles & 4 twin wheel castors	36.5kg

Customised Options

EndoTrans can be customised to suit your individual requirements using an array of optional accessory items.



Tray Kits

Standard Cleanascope® or Clear Cleanascope® Tray Kits can be added to allow storage of endoscopes or accessories (space allows up to 4 trays).

Product Ref	Product Description
106049S	Standard Cleanascope® Tray Kit includes 1 tray, 2 runners & fixings
106049C	Clear Cleanascope® Tray Kit includes 1 tray, 2 runners & fixings
104456	Clear Cleanascope® Tray Lid
104455	Clear Cleanascope® Tray
103615	Standard Cleanascope® Tray Lid
103614	Standard Cleanascope® Tray



Shelf Kit

Additional shelf kits can be added to provide further space for equipment units.

Product Ref	Product Description
106050	Shelf kit includes 1 x shelf & fixings

Accessory Items



Product Ref: 103339

Product Description:
LCD Monitor Post
Max weight: 20kg



Product Ref: 103342

Product Description:
Extension Post



Product Ref: 103922

Product Description:
Articulating Arm for LCD Monitor with Short Extension (500mm)
Max weight: 14.5kg



Product Ref: 104452

Product Description:
Twin Scope hanger



Product Ref: 103927

Product Description:
Articulating Arm for LCD Monitor with Long Extension (670mm)
Max weight: 14.5kg



Product Ref: 106052

Product Description:
Sterile Water Bottle Holder



Product Ref: 104454

Product Description:
co2 Gas Bottle holder (multi-size)



Product Ref: 101988

Product Description: Isolation Transformer Model EIT 2,000 VA



Product Ref: 104453

Product Description:
I.V Pole

+44 (0)1702 291878
+44 (0)1702 290013
info@medinn.co.uk
www.medicalinnovations.co.uk

Campfield Road, Shoeburyness, Essex, SS3 9BX, UK

