

Histolock™
resection device

Histo/ock
resection device



US endoscopy

listening...and delivering solutions®

Designed to grab and resect difficult flat lesions during endoscopic mucosal resection...


The Histolock™ resection device features a stiff monofilament wire to help recruit and resect flat tissue during complex polypectomy procedures.

The Histolock™ resection device features a...

twisted wire that lays flat on the mucosa during resection

torsion spring that helps maintain an open shape after multiple uses

mini snare size that offers easy maneuverability for controlled placement



For more information on the Histolock™ resection device, contact US Endoscopy today.

product number	description	catheter diameter (mm)	catheter length (cm)	snare size (width x length)	sterile	units/box
BX00711871	Histolock™ resection device	2.4	230	14mm x 27mm	yes	3
BX00711117	Histolock™ resection device	2.4	230	14mm x 27mm	yes	10

Distributed by:



listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com
761004 Rev. A

US Endoscopy is a subsidiary of STERIS Corporation.

© US Endoscopy or its affiliates. All rights reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

Follow us:    