


# INFLATION DEVICES

## High pressure syringes with manometer

- wide range of products
- wide pressure range
- ergonomic form
- clear drain, to facilitate observation of air bubbles



**BALTON**  
since  
1980

**BALTON**<sup>®</sup>

CE 0197

# INFLATION DEVICES

## High pressure syringes with manometer

Wide range of products for applications in:

- coronary angioplasty
- peripheral angioplasty
- valvuloplasty
- biliary ducts plasty
- gastrointestinal tract plasty

Reference:  
**SW 20**

pressure range: 0-30 atm
capacity: 20 cc
ergonomic form
easy and quick balloon deflation
90° rotating manometer
clear drain, to facilitate observation of air bubbles

Reference:  
**SW 25**

pressure range: 0-30 atm
capacity: 25 cc
ergonomic form
easy and quick balloon deflation
clear drain, to facilitate observation of air bubbles

Application: coronary angioplasty.

Reference:  
**SW 30**

pressure range: 0-30 atm
capacity: 30 cc
fluorescent indicator
easy and quick balloon deflation
clear drain, to facilitate observation of air bubbles

Application: coronary angioplasty, peripheral angioplasty, valvuloplasty.

Reference:  
**SW 60**

pressure range: 0-15 atm
capacity: 60 cc
ergonomic form
easy and quick balloon deflation
clear drain, to facilitate observation of air bubbles

Application: gastrointestinal tract plasty, valvuloplasty.

**Caution:**

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information on herein; these changes will be incorporated in the new editions of the publication. Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time. This publication is only a commercial information and cannot be handled the commercial offer in the meaning of the art.66 §1 of Polish Civil Code.

**BALTON**®

03-152 Warszawa, POLAND, ul. Modlińska 294  
tel.: +48 (22) 597 44 00, fax: +48 (22) 597 44 44  
e-mail: balton@balton.pl, www.balton.pl