


Lariat™
snare




Address various polyp sizes with one snare...


The multifunctional design of the Lariat™ snare helps address a variety of resection needs including piecemeal resection and multiple polyps.

The Lariat™ snare features...

- Three distinct sizes and shapes eliminating the need to exchange snares
- Snare loop design maintains its shape and integrity after multiple resections
- Stiffer snare wire minimizes distal deflection



Stage 1:
30mm oval



Stage 2:
10mm hex



Stage 3:
6mm diamond

For more information on the Lariat™ snare, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	approximate snare size (width x length)	units per box
00711119	Lariat™ snare	2.4	230	oval shape: 30x55mm hex shape: 10x28mm diamond shape: 6x15mm	10

Distributed by:



listening...and delivering solutions™

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

760931 Rev. C