

Double lumen tube distally tapered from 2.2 mm to 1.8 mm in treatment area.

double lumen



Tapered Tip

- for guide wire .035 inch



tapered tip

cutting length	20 mm	25 mm	30 mm
	OK 105 2220DL	OK 105 2225DL	OK 105 2230DL

double lumen



Standard Tip

- for guide wire .035 inch



standard tip

cutting length	20 mm	25 mm	30 mm
	OK 100 2220DL	OK 100 2225DL	OK 100 2230DL

double lumen



Precut-Type

- for guide wire .035 inch



precut tip

cutting length	20 mm	25 mm	30 mm
	OK 103 2220DL	OK 103 2225DL	OK 103 2230DL



Fr. 7



Fr. 5

134 °C
273 °F

Double lumen tube distally tapered from 2.2 mm to 1.8 mm in treatment area.

double lumen



Needle-Knife, Pointed-Type

- for guide wire .035 inch

stop at 6 mm

OK 110 18N-DL



Fr. 5

134 °C
273 °F

single lumen



Needle-Knife, Pointed-Type

- no guide wire access

stop at 6 mm

OK 110 18N



Needle-Knife, Blunt-Type

- no guide wire access

stop at 6 mm

OK 110 18SP



reusable

Needle-Knife Sphincterotomes

Fr. 7

134°C
273 F


single lumen



B-II-Type

- for guide wire .035 inch

OK 112 22FD



Fr. 5

134°C
273 F

single lumen



B-II-Type

- no guide wire access

OK 112 18

