








polypectomy snarees



An extensive line of polypectomy snares...

From cold snaring to advanced polypectomy, US Endoscopy offers a product solution to fit a variety of clinical needs.

product information

	product number	description	needle gauge/ projection	sheath diameter (mm)	length (cm)	approximate snare size (cm)	units/ box
	BX00711018	Short Throw™ snare - mini	-	2.4	230	1.25 x 3.0	10
	BX00711019	Short Throw™ snare - standard oval	-	2.4	230	2.5 x 5.4	10
	BX00711101	Short Throw™ snare - standard oval enteroscopy	-	2.4	350	2.5 x 5.4	5
	BX00711085	iSnare® system - oval	23 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711086	iSnare® system - oval	25 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711088	iSnare® system - hexagonal	25 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711089	iSnare® system - hexagonal	23 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711113	Traxtion® hexagonal snare	-	2.4	230	2.5 x 4.0	10
	BX00711115	Exacto® cold snare	-	2.4	230	9.0mm	10
	BX00711117	Histolock™ resection device	-	2.4	230	14mm x 27mm	10
	BX00711871	Histolock™ resection device	-	2.4	230	14mm x 27mm	3

Distributed by:

NEW!



	product number	description	needle gauge/ projection	sheath diameter (mm)	length (cm)	approximate snare size (cm)	units/ box
NEW!	BX00711872	Histolock™ resection kit - Histolock™ resection device Carr-Locke injection needle Mio® medical device organizer eTrap® polyp trap	-	-	-	-	3
NEW!	BX00711873	Histolock™ resection kit - Histolock™ resection device Carr-Locke injection needle Roth Net® retriever - polyp Mio® medical device organizer eTrap® polyp trap	-	-	-	-	3
NEW!	BX00711119	Lariat™ snare	-	2.4	230	oval: 30x50mm hex: 10x28mm diamond: 6x15mm	10
	BX00711120	Rotator™ snare - standard	-	2.4	230	2.5 x 5.4	10
	BX00711121	Rotator™ snare - mini	-	2.4	230	1.25 x 3.0	10

All active cord connections are Olympus style





US Endoscopy is a subsidiary of STERIS Corporation.
© US Endoscopy or its affiliates. All rights reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

5976 Heisley Road
Mentor, OH 44060 USA
phone +1 440 / 639-4494
fax +1 440 / 639-4495
email global@usendoscopy.com

www.usendoscopy.com

Follow us:    