



Combining the cleaning benefits of a bristled brush and squeegee into one device...

The DuoSwift™ combination squeegee brush features a tapered nylon brush on the leading end and a series of squeegees on the trailing end to effectively remove debris from flexible endoscopes while adhering to industry guidelines.



The DuoSwift™ combination squeegee brush features a...

- bristled brush to help loosen and remove bioburden with friction
- squeegee to pull liquid and remaining residue through the endoscope channel
- distinct design to help maximize cleaning touches within compatible endoscope channels

For more information on the DuoSwift™ combination squeegee brush, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	bristle diameter (mm)	squeegee diameter (mm)	channel size (mm)	units/box
BX00711618	DuoSwift™ combination squeegee brush	1.7	230	5-7	2.5-5.2	2.8-7.0	50
BX00711619	DuoSwift [™] combination squeegee brush with valve/control head cleaning brush	1.7	230	5-7 10 (valve/ conrol head)	2.5-5.2 (squeegee)	2.8-7.0	50

Distributed by:



listening...and delivering solutions®

5976 Heisley Road Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

