


Traxtion™
snare




*uniquely designed,
stiff .018" diameter
braided wire*



listening... and delivering solutions®

One of the stiffest braided wire snares...

Arguably one of the stiffest snares on the market today, the Traxtion™ snare is ideal for the resection of *flat* or *sessile lesions*. The *stiff braided wire* provides exceptional grip of the mucosa and retains its shape even after multiple resections.



Piecemeal resection of a sessile lesion using the Traxtion™ snare

Many snares tend to slip off of flat lesions and require repeated placement of the snare to complete the resection. With the *Traxtion™ snare*, you are able to maintain your grip of the mucosa without compromising your technique.

The Traxtion™ snare features...

- A uniquely designed, stiffer braided snare wire at .018" diameter requiring *no changes to electrocautery settings*
- The ability to retain its shape and integrity with multiple resections
- Memory points for better grip of mucosa

For more information on the Traxtion™ snare, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	approximate snare size (cm)	active cord style	unit per box
00711113	Traxtion™ snare	2.4	230	2.5 x 4.0	Olympus	10

Distributed by:



*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

760375 Rev B