



CapsoVision

empowering through innovation

See More. Care Better.

Enhance your practice with the
hassle-free CapsoCam® Plus System.

The CapsoCam® Plus System is a sophisticated, user-friendly capsule endoscopy solution that provides a **full 360° image** of the small bowel mucosa - the most comprehensive capsule endoscopy visualization available on the market!

The thorough imaging of the CapsoCam Plus System is enhanced by its wire-free design, which allows for flexibility in patient scheduling, while also providing unprecedented patient freedom and comfort during the examination.



CapsoView™




CapsoRetrieve®



CapsoAccess

CapsoAccess



Capso **CAM PLUS**

How it Works



CapsoRetrieve®

CapsoAccess®

CapsoView®



1 Patient swallows the small, vitamin-sized CapsoCam® Plus.

2 Patient retrieves the CapsoCam® Plus with the CapsoRetrieve® tool.

3 Patient data is rapidly downloaded with the CapsoAccess® system.

4 Physician quickly plays back and reviews patient data with the CapsoView® application.



Streamlined Procedure Saves Time and Money

CapsoCam® Plus captures a 360° panoramic lateral view of the small bowel using an on-board storage system that eliminates the need for belts, wires, leads and sticky pads.

Removing the need for external receiver equipment, the CapsoCam® Plus System comes with zero capital investment. Additionally, administration of a CapsoCam® Plus capsule takes less time and doesn't require patients to return bulky equipment that needs charging and maintenance, allowing for more efficient use of staff and exam room time.



Wire-free Design is Patient Preferred

CapsoCam® Plus patients are completely unburdened by uncomfortable external equipment. During the procedure, they can go about their daily activities while experiencing unprecedented comfort, convenience, and freedom.

The patient simply passes the capsule with a normal bowel movement, retrieves it using a specially designed retrieval kit, CapsoRetrieve®, and returns it to the clinic by mail or in person.



Fast and Easy Data Download is Hassle-Free

The CapsoAccess® capsule data access system allows you to seamlessly extract *in vivo* images and other data from the CapsoCam® Plus capsule after the procedure. To transfer data from the capsule to your computer, simply place the capsule into the CapsoAccess® system and close the lid.



User-friendly Software Makes Reading Faster and Easier

CapsoView® software provides a convenient, easy-to-use diagnostic user interface for viewing the CapsoCam® Plus images – with several playback options and simple, efficient report generation.



To learn more or to set up a demonstration, contact your local CapsoVision representative.

CapsoVision, Inc. ■ 18805 Cox Ave # 250, Saratoga, CA 95070
p. 408 624 1488 ■ f. 408 370 4795 ■ info@capsovision.com ■ www.capsovision.com

© 2016, CapsoVision, Inc., 01/2016 PN ADV-2554 Rev. B OUS